



Hoists & Accessories

Moving and handling solutions



The name SLK stands for quality
“Made in Germany”.
Service and customer satisfaction

SLK is a German developer and manufacturer of medical products and rehabilitation appliances.

With the SLK brand now as part of the Wellell global family, our customers can enjoy an even wider range of quality & innovative medical solutions enabling people around the world to enjoy a healthy and carefree life.

BUILT FOR PATIENT & CAREGIVER

All hoists are designed and built with both patient and caregiver in mind, perfected in co-operation with physicians, nursing staff and our partners.

The hoists are suitable for hospitals, care facilities and home care. The systems feature high quality and extra durable materials.

User Specification

	SLK Carry 185 Classic	SLK Multy	SLK Carry Compact	SLK Eazy-up
Max User Weight	185 kg/ 410Lbs/ 29 St	185 kg/ 410Lbs/ 29 St	160 kg/ 350Lbs/ 25 St	160 kg/ 350Lbs/ 25 St
User Height	210 cm/ 2.1 m / 6ft 10"	210 cm/ 2.1 m / 6ft 10"	170 cm/ 1.7 m / 5ft 7"	170 cm/ 1.7 m / 5ft 7"



UNIVERSAL HOISTS



ACTIVE HOISTS



ACCESSORIES & SPARES



Be well, Live well
Wellell

SLK Carry 185 Classic

Ideal solution for home care



Key Features & Benefits

- ✓ Soft-start/soft-stop
- ✓ Gentle pick up and lowering of the patient
- ✓ CPR Function (as standard)
- ✓ Battery indicator light on handset
- ✓ Comfort-control with display and readout function
- ✓ Raising & lowering via control box possible
- ✓ Lifting from the floor possible
- ✓ High battery power
- ✓ Highly manoeuvrable in small areas
- ✓ Hook mount with spring safety device
- ✓ 360-degree rotation of the spreader bar in highest possible position
- ✓ Minimises moving & handling risk to the caregiver
- ✓ Foldable without tools/ transportable
- ✓ Very quiet operation
- ✓ Standard mast and chassis are removable without tools



Standard

- SLK 2-Point Spreader Bar

Optional

- SLK 4-Point Spreader Bar X: 8060
- SLK 4 Point Spreader Bar H: 8064
- Electric Spreading : 8003



SLK 2-Point Spreader Bar




SLK 4-Point Spreader Bar X



SLK 4 Point Spreader Bar H

Specification (Hoist Product Code: 8002)

Manual Spreading	
Length	118 cm (46.4")
Height	137 cm (54")
Width	61 cm (24")
Leg Height	11 cm (4.3")
Lift Height	81-203 cm (31.8"-80")
Rotating Radius	132 cm (52")
Spread	51-94 cm (20"-37")
Total Weight	37 kg (81.5Lbs)
Castor Size	7.5 cm - 10 cm
Operating Noise	55 dB(A)



Max User Weight	185 kg/ 29 St/ 410Lbs
Max User Height	210cm/ 2.1m/ 6ft 10"
Material	<ul style="list-style-type: none"> • Powder coated steel • Ball-bearing swivel castors
Warranty	<ul style="list-style-type: none"> • 3 Years General • 6 Months Battery



www.wellell.com

SLK Multy

Active Hoist convertible to Universal Hoist



Wellell Group



MADE IN GERMANY

Key Features & Benefits

- ✓ Soft-start / Soft-stop
- ✓ Gentle pick up and settle of the patient
- ✓ CPR function (as standard)
- ✓ Telescopic stand mast
- ✓ Battery indicator light
- ✓ Comfort-control with display and readout function
- ✓ Handling via control box possible
- ✓ Lifting from the floor possible with Active Hoist Kit
- ✓ High battery power
- ✓ High manoeuvrability in a small area
- ✓ Hook mount with spring safety device (spreader bar only)
- ✓ 360-degree rotation of Spreader Bar in highest possible position (in universal version only)
- ✓ Protects the back of the caregiver
- ✓ Very quiet operation
- ✓ Foldable without tools / transportable
- ✓ Standard mast and chassis are removable without tools



Standard

- Footboard
- Shin Pad
- Active Lifting Arm

Optional


- SLK Multy Conversion Kit (Lifting Arm & Spreader Bar): 8010
- Electric Spreading : 8013



Specification - Active Hoist Configuration

(Hoist Product Code: **8008**)

Manual Spreading	
Length	118 cm (46.4")
Height	105 cm (41.3")
Width	61 cm (24")
Leg Height	11 cm (4.3")
Lift Height	85-165 cm (33.4"-65")
Rotating Radius	132 cm (52")
Spread	51-94 cm (20"-37")
Weight	40kg (88lbs)



Max User Weight	185 kg / 29 St / 410 Lbs
Max User Height	210cm / 2.1m / 6ft 10"
Material	<ul style="list-style-type: none"> • Powder coated steel • Ball-bearing swivel castors
Warranty	<ul style="list-style-type: none"> • 3 Years General • 6 Months Battery



www.wellell.com

SLK Multy Conversion Kit


for SLK Multy Active Hoist



Wellell Group



MADE IN GERMANY



Conversion Kit Standard

- Lifting Arm
- SLK 2-Point Spreader Bar

Optional

- SLK 4-Point Spreader Bar X: 8060
- SLK 4 Point Spreader Bar H: 8064
- Electric Spreader



SLK 2-Point Spreader Bar

SLK 4 Point Spreader Bar H

SLK 4-Point Spreader Bar X

Specification - with Universal Hoist Kit

(Hoist Product Code: **8010**)

Manual Spreading	
Length	118 cm (46.4")
Height	137 cm (54")
Width	61 cm (24")
Leg Height	11 cm (4.3")
Lift Height	81-203 cm (31.8"-80")
Rotating Radius	132 cm (52")
Spread	51-94 cm (20"-37")
Weight	37 kg (81.57Lbs)



SLK Multy Conversion Kit Universal

Max User Weight	185 kg / 29 St / 410 Lbs
Max User Height	210cm / 2.1m / 6ft 10"

SLK Carry Classic Compact

Universal Hoist with short chassis for narrow living spaces



Key Features & Benefits

- ✓ For children and patients up to 1.7m / 5ft 7"
- ✓ Particularly low front under-ride height of 83mm
- ✓ Soft-start/ Soft-stop
- ✓ Gentle pick up and settle of the patient
- ✓ CPR function (as standard)
- ✓ Battery indicator light on handset
- ✓ Comfort-control with display and readout function
- ✓ Handling via control box possible
- ✓ Lifting from floor level possible
- ✓ Easy to use in narrow spaces and doors
- ✓ Adjustable chassis spreading
- ✓ High manoeuvrability in a small area
- ✓ High battery power/ Li-Ion battery
- ✓ Hook mount with spring safety device
- ✓ 360-degree rotation of the Spreader Bar in highest possible position
- ✓ Protects the back of the caregiver
- ✓ Foldable without tools/ transportable
- ✓ Standard mast and chassis are removable without tools

MADE IN GERMANY



Standard

- SLK 2-Point Spreader Bar Compact H

Optional

- SLK 4-Point Spreader Bar X: 8060
- SLK 4 Point Spreader Bar H: 8064



SLK 2-Point Spreader Bar Compact H



SLK 4-Point Spreader Bar X



SLK 4 Point Spreader Bar H

Specification

(Hoist Product Code: **8014**)

Manual Spreading	
Length	103 cm (40.5")
Height	124 cm (48.8")
Width	55 cm (21.6")
Leg Height	8.3 cm (3.2")
Lift Height	71-193 cm (28"-76")
Rotating Radius	114 cm (45")
Spread	45-78cm (17.7"-30.7")
Weight	33 kg (72.7Lbs)



Max User Weight	160 kg/ 350 Lbs/ 25 St
Max User Height	170 cm, 1.7m, 5ft7" (67")

Material	<ul style="list-style-type: none"> • Powder coated steel • Ball-bearing swivel castors
Warranty	<ul style="list-style-type: none"> • 3 Years General • 6 Months Battery



www.wellell.com

SLK Eazy-up

The Active Hoist with fixed, particularly low chassis for driving under the seat



Key Features & Benefits

- ✓ Soft-start/ Soft-stop
- ✓ Gentle pick up and lower of the patient
- ✓ CPR function (as standard)
- ✓ Battery indicator light on handset
- ✓ Comfort control with display and readout function
- ✓ Easy to use in narrow spaces and doors
- ✓ High manoeuvrability in small area
- ✓ Easier standing up due to particularly deep footboard
- ✓ Adapted to the natural get up movement
- ✓ Safe transfer even in standing position
- ✓ Therapy-supporting, for preservation or increase mobility
- ✓ Protects the back of the caregiver
- ✓ Pluggable mast
- ✓ Standard mast and chassis are removable without tools

MADE IN GERMANY



Standard

- SLK Shin Pad

Optional

- SLK Narrow Calf Band: 8094
- Manual Spreading Chassis: 8016


Appropriate Slings

- SLK Active Sling Eco (size M/L)
- **SLK Active Sling (size S - XXL)**
- **SLK Special Active Sling (size S - L)**
- **SLK Special Active Sling with Seat Swing (size S-XXL)**

Specification

(Hoist Product Code: **8017** fixed chassis)

Manual Spreading	
Length	99 cm (39")
Height	114 cm (45")
Width	58 cm (22.8")
Leg Height	6.6 cm (2.5")
Lift Height	90-159 cm (35.4"-62.5")
Rotating Radius	90 cm (35.4")
Weight	34 kg (75 Lbs)




Max User Weight	160 kg/ 350 Lbs/ 25 St
Max User Height	170 cm, 1.7m, 5ft7" (67")

Material	<ul style="list-style-type: none"> • Powder coated steel • Ball-bearing swivel castors
Warranty	<ul style="list-style-type: none"> • 3 Years General • 6 Months Battery



www.wellell.com



MADE
IN
GERMANY



www.wellell.com



www.wellell.com

	SLK 2-Point Spreader Bar	SLK 4-Point Spreader Bar X	SLK 4-Point Spreader Bar H	SLK 4-Point Spreader Bar H (Compact)
Width	59 cm (23.2")	60 cm (23.6")	57 cm (22.4")	52 cm (20.5")

	SLK Carry 185 Classic	SLK Multy Universal	SLK Carry Compact	SLK Eazy-up
SLK 2-Point Spreader Bar 8053 (telegrey)	✓	-	-	-
SLK 2-Point Spreader Bar 8054 (silver)	✓	✓	-	-
SLK 2-Point Spreader Bar 8689 (telegrey)	-	-	✓	-
SLK 4-Point Spreader Bar X 8060 (silver)	✓	✓	✓	-
SLK 4-Point Spreader Bar H 8064 (silver)	✓	✓	✓	-
SLK 4-Point Spreader Bar H 8095 (telegrey)	-	-	✓	-

SLK Velcro Headrest



SLK Calf Band




SLK Adjustable Headrest Reinforced



SLK Toilet Sling



SLK Quick Transfer and Toilet Sling



SLK Special Active Sling



SLK Special Active Sling with Seat Swing



SLK Active Sling



SLK Seat Swing




SLK Universal Sling



SLK Universal Sling with Headset




SLK Back-Fix Sling




SLK Maxi Sling




SLK Hammock Seat Pan Sling



SLK Hammock Universal Sling



SLK Comfort Sling



SLK 2-Point Spreader Bar

The SLK 2-Point Spreader Bar allows a particularly comfortable transport feeling. Slings cannot be unintentionally released due to its spring lock

- ID 8053 (telegrey) : for SLK Carry 185 Classic
- ID 8054 (silver) : for SLK Carry 185 Pro, Pro L, SLK Multy Universal and electrically spreadable models.
- ID 8689 (telegrey): for SLK Carry Compact Universal Hoist



SLK 4-Point Spreader Bar X

The SLK 4-Point Spreader Bar X allows a particularly comfortable transport feeling. Slings cannot be unintentionally released due to its spring lock

- ID 8060 (silver) : for SLK Carry Product family, SLK Multy Universal and the electrically spreadable models



SLK 4-Point Spreader Bar H

The SLK 4-Point Spreader Bar H allows a particularly comfortable transport feeling. It is especially suitable for plus size patients. Slings cannot be unintentionally released due to its spring lock

We recommend the use of the SLK 4-Point Spreader Bar H inter alia to:




























- Plus size patients
- Sensitive patients
- Patients with stretch spasticity issues










- ID 8064 (silver) : for SLK Carry Product family, SLK Multy Universal and the electrically spreadable models

- ID 8064 (telegrey) : for SLK Carry Compact Universal Hoist









Hoist Components

Article IDS of Spare Parts		SLK Multy				Eazy Up		Carry				
Spare Part *1	Drawing	Active	Universal	Conversion kit Active	Conversion Kit Universal	Eazy-up Fix	Eazy-up flex	185 eco	185 Classic	185 Pro	185 Pro L	Compact
User Guide (EN)		8293	8293	8293	8293	8292	8292	8304	8304	8304	8304	8304
Foot, left		8656	8656	-	-	-	8758	8654	8654	8656	8658	8653
Foot right		8662	8662	-	-	-	8759	8661	8661	8662	8663	8660
Central Tube		8736	8736	-	-	-	8760	8665	8665	8666	8666	8664
Mast		-	-	-	-	8761	8761	8669	8670	8672	8674	8668
Pusher Bar		8678	8678	-	-	8662	8662	8677	8678	8678	8678	8678
Lifting Arm		-	-	-	-	8763	8763	8680	8681	8682	8683	8679
2-Point Spreader Bar		-	8054	-	8054	-	-	8053	8053	8054	8054	8689
Pedal		8652	8652	-	-	-	8650	8651	8651	8652	8652	8650
Cinematics		-	-	-	-	-	-	-	-	-	-	-
Spreader Bar		8692	8692	-	-	-	8691	8692	8692	8692	8692	8691
Caster front		8696	8696	-	-	8764	8693	8694	8695	8696	8696	8693
Caster rear		8700	8700	-	-	8697	8697	8698	8699	8700	8700	8697
Service Set*2		-	-	-	-	-	-	-	-	-	-	-
Inner Mast		8723	8723	-	-	-	-	-	-	-	-	-
Outer Mast		8725	8725	-	-	-	-	-	-	-	-	-
Gas Pressure-spring		8731	8731	-	-	-	-	-	-	-	-	-
Lifting Arm Base		8721	8721	-	-	8763	8763	-	-	-	-	-
Lifting Arm Active		8720	-	8720	-	8757	8757	-	-	-	-	-
Lifting Arm Universal		-	8738	-	8738	-	-	-	-	-	-	-
Service Set *3		8705	8705	-	-	8765	8765	8704	8704	8704	8704	8704
Footboard flat		8728	-	8728	-	-	-	-	-	-	-	-
Footboard deep		8729	-	8729	-	-	-	-	-	-	-	-
Shin Support metal		8733	-	8733	-	8733	8733	-	-	-	-	-
Shin Support Soft		8735	-	8735	-	8735	8735	-	-	-	-	-
Shin Support Pad		8734	-	8734	-	-	-	-	-	-	-	-
Base		-	-	-	-	8755	-	-	-	-	-	-

Drive Components Type I*3

Article IDS of Spare Parts		SLK Multy				Eazy Up		Carry				
Spare Part *1	Drawing	Active	Universal	Conversion kit Active	Conversion Kit Universal	Eazy-up Fix	Eazy-up flex	185 eco	185 Classic	185 Pro	185 Pro L	Compact
Hand-control		8636	8636	-	-	8636	8636	8636	8636	8636	8636	8636
Charger		8633	8633	-	-	8633	8633	8633	8633	8633	8633	8633
Battery Pack		8632	8632	-	-	8632	8632	8632	8632	8632	8632	8632
Control Box		8631	8631	-	-	8631	8631	8631	8631	8631	8631	8631
Actuator		8637	8637	-	-	8637	8637	8637	8637	8637	8637	8637
Base Expander Motor		-	-	-	-	-	-	-	-	-	-	-
Safety Box		8640	8640	-	-	-	-	-	-	-	-	-

Drive Components Type II*3

Article IDS of Spare Parts		SLK Multy				Eazy Up		Carry				
Spare Part *1	Drawing	Active	Universal	Conversion kit Active	Conversion Kit Universal	Eazy-up Fix	Eazy-up flex	185 eco	185 Classic	185 Pro	185 Pro L	Compact
Hand-control		8445	8445	-	-	8445	8455	8446	8445	8445	8445	8445
Hand-control Magnet		8455	8455	-	-	8455	8455	-	8455	8455	8455	8455
Charger		8450	8450	-	-	8450	8450	8450	8450	8450	8450	8450
Battery Pack		8441	8441	-	-	8441	8441	8442	8442	8441	8441	8442
Control Box		8444	8444	-	-	8444	8444	8443	8444	8444	8444	8444
Safety Box		8722	8722	-	-	-	-	-	-	-	-	-
Actuator		8396	8396	-	-	8396	8396	8449	8396	8396	8396	8396
Hand-control Clip		-	-	-	-	-	-	8456	-	-	-	-

*1 Each spare part is delivered with the small parts that are required for assembly.
 *2 The service set includes all screws, nuts, stoppers and other parts.
 *3 Information about the drive system is given on the relevant type labels

HEADQUARTERS

Wellell Inc.

No. 9, Minsheng St., Tucheng Dist.,
New Taipei City 23679, Taiwan
T +886 2 2268 5568
F +886 2 2268 6525
marketing@wellell.com

PAN-ASIA

CHINA

Wellell (Kunshan)Co., Ltd.
No. 1368, ZiZhu Road, Kunshan,
Jiang Su Province, China
T +86 512 8617 1660 Ext 1150
F +86 512 8617 8663
customer.cn@wellell.com

THAILAND

Wellell (Thailand)
No.111/ 152 Soi Pho Kao Yak 19
Nawamin, Beung Kum,
Bangkok 10240,
Thailand
T +662 2018 2140
sales.thai@wellell.com

PAN-EUROPE

SPAIN

Apex Medical S.L.
Elcano 9, 6ª planta 48008
Bilbao. Vizcaya. Spain
T +34 94 4706408
F +34 94 4706409
info@apexmedical.es

UK

Wellell UK Limited
Unit 33, Great Western Business Park
McKenzie Way, Worcester, WR4 9GN
United Kingdom
T +44 (0) 1905 774 695
sales.uk@wellell.com

FRANCE

Wellell France S.A.S
4 Boulevard de la Chanterie
49124 Saint-Barthelemy d' Anjou
France
T +33 2 72 79 27 13
F +33 2 85 73 07 67
info.france@wellell.com

GERMANY

SLK Vertriebsgesellschaft GmbH
Am Herdicksbach 18, 45731 Waltrop
Germany
T +49 (0) 2309 91545 0
F +49 (0) 2309 91545 999
info@slk-gmbh.de
www.slk-gmbh.de

AMERICAS

USA

Wellell America Corp.
927 Mariner Street,
Brea, CA 92821, USA
T +1 714 671 3818
F +1 714 494 8173
sales.usa@wellell.com