

## Icon Active

Delivering ease of care and comfort



## Icon Active Delivering ease of care and comfort.

Designed using quality materials, this electronically controlled community bed offers safety, support and comfort.

This stylish bed benefits from the additional design features of an adjustable backrest, leg-lift and overall bed height adjustable from 24cm to 64cm offering comfort, flexibility and ease of transfer.

Ideal for non-acute settings where patient comfort and mobility is paramount, this community bed is designed to support and assist not just the patient but the caregiver where needed.



## Key Features

- Easy assembly and installation
- Electric Wooden bed with adjustable head and foot sections
- Height adjustable from 24cm to 64cm
- Auto profiling feature
- Trendelenburg/Reverse Trendelenburg feature
- Electric height adjustment
- Electric knee break
- Easy clean and maintenance
- 3 year warranty on metalwork and actuators

## Specifications

Maximum user weight	:180KG
Platform height	: 24cm-64cm
Mattress base width	: 89 cm- length 202cm
Foot and head panel height from mattress base	: 47cm
Castor Diameter	: 32 (PVC) WH001
Backrest angle	: 30° degrees
Trendelenburg/Reverse Trendelenburg	
Extendable option	
Lockable handset	
Static foam and dynamic mattresses available	

UK

Apex Medical Limited Unit 33, Great Western Business Park McKenzie Way, Worcester, WR4 9GN T: +44 (0)1905 774 695 sales@apexmedicalcorp.co.uk www.apexmedicalcorp.co.uk



www.apexmedicalcorp.co.uk

© Apex Medical Corp. 2019. All rights reserved. Specifications may change without notice.

Distributed by