

Icon Floor & Icon Floor Plus

Minimising the risk of fall and potential harm



Icon Floor & Icon Floor Plus

Minimising the risk of fall and potential harm

The electronically controlled Icon Floor bed range is designed specifically with the occupant in mind, combining quality construction and functionality to assist with care, as well as minimising the risk to the occupant of falling and causing injury by using the low height function.

The stylish bed has additional features including backrest, leg-lift and an overall bed height adjustable from 10cm to 70cm to further aid mobility and care.

Using a combination of these functions, the Icon Floor range can also be placed into a comfortable profiling position for daily activities such as watching television or interacting with visitors making it the perfect profiling bed for multiple care environments including home, residential and nursing care.



Key Features

- Minimum mattress platform height: 10cm from the floor
- Easy to use handset
- 2 stage lowering to prevent entrapment
- Electric angle adjustment of the foot section
- Auto contour function
- Auto regression function
- Easy to clean and maintain
- Aesthetically designed head/foot boards for enhanced visibility and a fully enclosed foot board option
- Independently braked castors
- Quick and easy assembly
- 3 year warranty on metalwork and electrical actuators



Wellell UK Limited Unit 33, Great Western Business Park McKenzie Way, Worcester, WR4 9GN T: +44 (0)1905 774 695 sales.uk@wellell.com



www.wellell.com

Specifications

Maximum user weight	:180KG
Lowest platform height	: 10cm from floor
Maximum height position	: 70cm
Mattress platform	:89cm-200cm
Mattress platform extension	: 15cm (total = 215cm)
Castor diameter	: 3" (PVC) WH001
Backrest angle	: 73°
Calf section angle (when thigh section raised)	: 25°
Thigh section raised angle	: 30°
Auto regression function	
Trendelenburg/reverse trendelenburg	
Auto profile function	
Anti-entrapment function when lowering below 30 cm	

Static foam and dynamic mattresses available

Distributed by

Icon Floor Plus