



Vac Plus

Portable Suction Unit

The Vac series pump is APEX Medical's most recent generation of portable suction units and has been designed to combine function with aesthetic appeal and meet the needs of today's hygiene and infection control requirements. The Vac Plus unit is lightweight, compact and can be used in either AC or DC mode. This unit is available with a rechargeable battery, remote controller, and car adaptor.

Vac Plus

Portable Suction Unit

Features ▼

- Smooth hygienic finish, water resistant, easy to clean and maintain.
 - For home and outdoor use.
 - Optional accessories available
- Single-patient use canister available in 800ml or 1800ml.
- Filter and Tubing.



Integrated Canister Holder



Portable design with ergonomic easy-to-carry handle.



Pressure regulator and easy-to-read analogue pressure gauge.



Optional attached remote controller for additional operational convenience.



Single-patient use canister available in 800ml or 1800ml.

Specification

Electrical Rating	AC 100 - 240V, 50/60 Hz ; DC 12V
Power Consumption	40W
Maximum Pressure	600 mmHg (adjustable)
Maximum Flow Rate	24 LPM
Net Weight	3.3 kg
Dimension	352 (L) x 206 (W) x 192 (H) mm
Sound Level	< 55 dB



Distributed by

Apex Medical Corp.
No. 9, Min Sheng St., Tu-Cheng,
New Taipei City, 23679, Taiwan
T +886 2 2268 3100
F +886 2 2268 6525
marketing@apexmedicalcorp.com
www.apexmedicalcorp.com
Print-2017/01 All rights reserved
Specifications may change without notice.