

HEADQUARTERS

Apex Medical Corp.  
No. 9, Min Sheng St., Tu-Cheng,  
New Taipei City, 23679, Taiwan  
T +886 2 2268 5568  
F +886 2 2268 6525  
marketing@apexmedicalcorp.com

PAN-ASIA

**TAIWAN**  
Apex Medical Corp.  
No. 9, Min Sheng St., Tu-Cheng,  
New Taipei City, 23679, Taiwan  
T +886 2 2268 5600 /0800 054 133  
F +886 2 2268 6526  
customer@apexmedicalcorp.com

**STURDY INDUSTRIAL**  
168, Sec.1, Zhongxing Road, Wugu District,  
New Taipei City 24872, Taiwan  
T +886 2 2981 5245  
F +886 2 2984 8208  
foreign@sturdy.com.tw  
www.sturdy.com.tw

**CHINA**  
Apex Medical (Kunshan) Corp.  
No. 1368, ZiZhu Road, Kunshan,  
Jiang Su Province, China  
T +86 512 8617 1660 Ext 1150  
F +86 512 8617 8663  
customer.cn@apexmedicalcorp.com

**INDIA**  
Apex Medicalcorp India Pvt. Ltd.  
303, Third Floor Block 4A DLF  
Corporate Park Phase-III,  
Gurgaon, Haryana  
T +91 11 6654 4138  
F +91 11 6654 4052  
enquiries.india@apexmedicalcorp.com

**THAILAND**  
Apex Medical Thailand  
65/100 Chamnanphenjati Business Center Bldg.  
11/F Rama 9 Road, Huaykwang  
Bangkok 10320  
T +662 948 4848  
F +662 948 4170  
sales.thai@apexmedicalcorp.com

**MIDDLE EAST • AFRICA**  
29, Misr Liltameer Building,  
Stree#10, AREA#8, Sheraton Building  
Heliopolis Cairo, Egypt.  
T +202 2 2671303  
F +202 2 2678649  
raafat.wahby@apexmedicalcorp.com

PAN-EUROPE

**SPAIN**  
Apex Medical S.L.  
Elcano 9, 6ª planta 48008  
Bilbao, Vizcaya, Spain  
T +34 94 4706408  
F +34 94 4706409  
info@apexmedical.es

**UK**  
Apex Medical Limited  
Unit 33, Great Western Business Park  
McKenzie Way, Worcester ,WR4 9GN  
United Kingdom  
T +44 (0) 190 5774 695  
sales@apexmedicalcorp.co.uk

**FRANCE**  
Apex Medical France  
4 Boulevard de la Chanterie  
49124 Saint-Barthelemy d'Anjou  
France  
T +33 2 72 79 27 13  
F +33 2 85 73 07 67  
info.france@apexmedicalcorp.com

**GERMANY**  
SLK Vertriebsgesellschaft GmbH  
Lindenhorster Straße 38-40 44147  
Dortmund  
T +49 (0)231-925360-175  
F +49 (0)231-925360-29  
info@slk-gmbh.de  
www.slk-gmbh.de

AMERICAS

**USA**  
Apex Medical USA  
927 Mariner Street,  
Brea, CA 92821, USA  
T +1 714 671 3818  
F +1 714 494 8173  
sales.usa@apexmedicalcorp.com



ACS Turn  
Low Air Loss Turning Air Mattress

- Suitable for patients in high to very high risk of pressure injuries
- Up to 40° continuous lateral rotation to improve secretion clearance
- True low air loss therapy offers better skin microclimate management
- Designed for use in the acute care sector



Care for a Healthy Life



Distributed by



# ACS Turn

## Low Air Loss Turning Air Mattress



### Improving Pulmonary Function and Oxygenation with 40° Continuous Lateral Rotation


ACS Turn with true low air loss system can provide an up to 40 degree continuous lateral rotation which has been clinically approved to be effective in reducing the pulmonary complications.\*<sup>1</sup>

References:  
\*1. Wilson, D. et al., Wong's Clinical Manual of Paediatric Nursing-E-Book. Elsevier Health Sciences, 2014

Shift-Stop Technology  
Ensures rotational therapy  
safe and effective with  
special designed cells



Drainage Tubing  
Placement  
Allows drainage tubing  
to pass through easily



4-way Stretch Cover

- Antimicrobial
- Low friction and shear
- Vapour permeable, water resistant
- High-flame retardant



Water Resistance



Flame Resistance



Biocompatibility



Fungistatic



4-Way Stretch



Vapour Permeable



Silver-ion Treated



True Low Air Loss  
Better skin microclimate  
management to minimize  
tissue deformation and  
prevent skin breakdown



Visual and Audible Alarm

- Low Pressure
- Power Failure
- Tech Support



True Low Air Loss with  
Blower System



Cable  
Placement  
Hook



Care-Assist Foam  
Base Pocket  
Reserved 2" foam base  
pocket and 10" side  
bolster for additional  
patient safety




2-way Zipper  
For convenient  
removal of covers  
for cleaning and  
maintenance



Obvious CPR Connector

- CPR < 15 seconds
- Storage vacuum mode



Continuous Lateral Turning

- Effortless to facilitate a turn
- Improving bronchial clearance for patients



Turning Options

- Selectable turning angles and dwell times can be tailored to meet patient's individual need



Max Firm Mode

- For quick patient transfer or nursing intervention

### Specification

| Model    | ACS Turn                  | 10" Replacement     |
|----------|---------------------------|---------------------|
| Pump     | Dimensions                | 15.8 x 38 x 26.8 cm |
|          | Weight                    | 6 kg / 13.2 lbs     |
|          | Case Material             | Flame retardant ABS |
|          | Supply Voltage            | AC 220~240, 50Hz    |
| Mattress | Dwell Time                | 5, 30, 60 min       |
|          | Dimensions                | 200 x 90 x 25.4 cm  |
|          | Cell Height               | 8" x 20 Cells       |
|          | Weight                    | 10.1 kg / 22.2 lbs  |
|          | Cover Material            | 4-way stretch       |
|          | Cell Material             | Nylon/ PU           |
|          | Maximum Patient Weight    | 180 kg / 397 lbs    |
|          | Flame Retardant Standards | EN 597-1 & EN 597-2 |

Specifications are subject to change for improvement without notice