

## Sedens 410 Pressure Relieving Seat Cushion

Pressure Relief for a 24 Hour Care Package

The Sedens 410 is a light weight and portable seat cushion designed specifically to be used in chairs along with APEX dynamic mattress system pumps. It's able to extend pressure relief therapy to sit position that provides patient 24 hour pressure area care, allowing caregivers to make versatile care plan for patient to prevent pressure injuries and increase patient's mobility.



# Sedens 410

Pressure Relieving Seat Cushion

### Providing 24 Hour Pressure Area Care for Patient

- Alternating cushion works with APEX dynamic mattress system pumps.
- Hybrid design with alternating therapy to relieve pressure from patient's buttock whilst reducing the risk of bottoming-out.



encased in to reduce risk of bottoming out.

Non-slip base with straps prevents the cushion from sliding out of position, giving peace of mind to the user.

#### Specifications

Cushion	Cushion Dimension	44 x 41 x 9 cm
	Cell Height	3″
	Weight	1.11 kg
	Cover Material	Stretch PU Polyester
	Cell Material	Nylon PU
	Base Material	Nylon PU
	Fire Retardant Rating	EN 597-1 & EN 597-2
	Maximum Patient Weight	No weight limit, yet the cushion must be properly sized to the individual

Specifications are subject to change for improvement without notice

#### **HEADQUARTERS**

Apex Medical Corp. No. 9, Min Sheng St., Tu-Cheng, New Taipei City, 23679, Taiwan T +886 2 2268 5568 F +886 2 2268 6525 marketing@apexmedicalcorp.com

#### UK

Apex Medical Limited Unit 33, Great Western Business Park McKenzie Way, Worcester, WR4 9GN . United Kingdom T +44 (0) 190 5774 695 sales@apexmedicalcorp.co.uk