

# Sedens 500

## Dynamic Seat Cushion

### Lightweight seat cushion with long-lasting battery power

The Sedens 500 aligns with the 2019 clinical guideline, featuring advanced cushions and back supports designed to prevent pressure injury and manage pressure effectively. Lightweight and portable, it's perfect for chairs or wheelchairs, offering exceptional postural stability with a solid foam base and medial thigh support. Its long-lasting lithium-ion battery provides up to 12 hours of uninterrupted use, enhancing convenience. Visual and audible alarms prioritize user safety, making the Sedens 500 a comprehensive and innovative solution for mobility challenges while promoting comfort and well-being.



#### Optimal Pressure Redistribution

Effective alternating mode and static therapy mode to fulfill individual needs, and it also offers 5 comfort levels and 3 cycle times option for selection.



#### Extended Outdoor Use

The device enables extended outdoor use with its 12-hour lithium-ion battery on a full charge.

# Sedens 500

Advanced mobility seat cushion: prolonged power supply and safety features



Specifications	Sedens 500			
<b>Pump</b>	Dimension	25 x 12.5 x 5.7 cm (9.8 x 4.7 x 2.2 in)		
	Weight	1 kg (2.2 lb)		
	Supply voltage	100-240V /50-60Hz		
	Therapy mode	Alternating/ Static		
<b>Cushion</b>	Dimension	Length	Width	Height
		41-46cm(16.1-18.1 in)	42-46cm(16.5-18.1 in)	10cm(4 in)
	Type	contour foam base with 4 air cells		
	Weight	1-1.2 kg (2.2-2.6 lb)		
	Top cover material	Poly/PU		
	Cell material	Nylon/PU		
	Base material	Nylon/PU		
Maximum patient weight	120kg (264 lb)			

Pump: water resistant standards (IP22); Mattress: flame retardant standards (EN597-1, EN597-2), RoHS, WEEE

